Inspired by
temperature

Externe Kùhlschlange
External cooling coil

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Externe Kühlschlange
External cooling coil

External cooling coil

This installation guide is a translation of the original installation guide.

VALID IN CONJUNCTION WITH:

Huber temperature control unit

Table of contents

V1.0.0en/21.01.20

1	Introduction	10
1.1	Safety	10
1.1.1	Symbols used for Safety Instructions	10
1.1.2	Representation of safety identifiers	11
1.1.3	Safety during commissioning	11
1.1.4	Extension of specified normal operation	11
1.2	Responsible bodies and operators – Obligations and requirements	11
1.2.1	Obligations of the responsible body	11
1.2.1.1	Proper disposal of resources and consumables	12
1.2.2	Requirements for operators	12
1.2.3	Obligations of the operators	12
1.3	General information	13
1.3.1	Description of workstation	13
2	Commissioning	14
2.1	Unpacking	14
2.2	Ambient conditions	14
2.3	Installation conditions	15
2.4	Recommended temperature control and cooling water hoses	15
2.5	Wrench sizes and torques	16
2.6	Temperature control units with cooler	16
2.7	Connecting an externally closed/open application	18
3	Function description	19
3.1	Function description of the accessory	19
3.1.1	General functions.....	19
4	Setup mode	20
4.1	Filling and draining the accessory	20
5	Normal operation	21
5.1	Counter-cooling	21
5.1.1	Start counter-cooling	21
5.1.2	End counter-cooling	21
6	Service/maintenance	22
6.1	Maintenance	22
6.1.1	Function check and visual inspection	22
6.1.2	Replacing temperature control or coolant hoses	22
6.1.2.1	Replacing temperature control hoses	22
6.1.2.2	Replacing coolant hoses	23
6.2	Thermal fluid inspection, replacement and circuit cleaning	23
6.3	Cleaning the surfaces	23
6.4	Decontamination/repairs	23
7	Shutting down	24
7.1	Safety instructions and basic principles	24

7.2	Switch-off	24
7.3	Draining the cooling water	25
7.4	Draining the accessory	25
7.5	Separating the accessory from the temperature control unit	25
7.6	Packing	26
7.7	Shipping.....	26
7.8	Disposal	26
7.9	Contact data	26
7.9.1	Telephone number: Customer Support	26
7.9.2	Telephone number: Sales	27
7.9.3	Email address: Customer Support.....	27
7.10	Certificate of Compliance	27
8	Annex	28
8.1	Technical data	28

Foreword

Dear Customer,

Thank you for choosing an accessory from Peter Huber Kältemaschinenbau SE. You have made a good choice. Thank you for your confidence.

Please read this installation guide carefully before commissioning. Strictly follow all notes and safety instructions.

Follow the installation guide with regard to transport, start-up, operation, maintenance, repair, storage and disposal of the accessory.

We fully warrant the accessory for the specified normal operation.

The component listed on page 5 is referred to in this installation guide as accessory and Peter Huber Kältemaschinenbau SE as Huber company or Huber.

Liability for errors and misprints excluded.

The following trademarks and the Huber logo are registered trademarks of Peter Huber Kältemaschinenbau SE in Germany and/or other countries worldwide: BFT®, CC®, Chili®, Com.G@te®, Compatible Control®, CoolNet®, DC®, E-grade®, Grande Fleur®, KISS®, Minichiller®, Ministat®, MP®, MPC®, Peter Huber Minichiller®, Petite Fleur®, Pilot ONE®, RotaCool®, Rotostat®, SpyControl®, SpyLight®, Tango®, TC®, UC®, Unical®, Unichiller®, Unipump®, Unistat®, Unistat-Pilot®, Unistat Tango®, Variostat®. The following trademarks are registered in Germany to DWS Synthesetechnik: DW-Therm®, DW-Therm HT®

1 Introduction

1.1 Safety

1.1.1 Symbols used for Safety Instructions

Safety instructions are marked by the below combinations of pictograms and signal words. The signal word describes the classification of the residual risk when disregarding the installation guide.

Denotes an immediate hazardous situation that will result in death or serious injuries.

Denotes a general hazardous situation that may result in death or serious injuries.

Denotes a hazardous situation that can result in injury.


Denotes a situation that can result in property material damage.

Denotes important notes and usable hints.

Safety information and procedure

The safety information in the installation guide is designed to protect the responsible body, operator and the equipment from damage. First inform yourself about any residual risks due to misuse before you start an operation.

1.1.2 Representation of safety identifiers

The following pictograms are used as safety identifiers. The table gives an overview of the safety identifier used here.

Identifier	Description
Mandatory sign	
	- Observe the instructions
Warning sign	
	- General warning sign - observe the instructions
	- Warning of electrical voltage
	- Warning of hot surface
	- Warning of flammable substances

1.1.3 Safety during commissioning

The following chapters are only relevant for accessory in connection with a HUBER temperature control unit and apply in addition to the operation manual of the temperature control unit used. If you have any questions regarding the installation guide, please contact our customer support (see page 26 in the section »Contact data«). The installation guide is to be kept for future use.

1.1.4 Extension of specified normal operation

The accessory is suitable for counter-cooling when properly installed at the temperature control unit. The accessory itself can **not** be used without being connected to the temperature control unit. Otherwise the intended use as described in the temperature control unit's operation manual applies.

1.2 Responsible bodies and operators – Obligations and requirements

1.2.1 Obligations of the responsible body

Store the installation guide near the accessory where it is easy to access. Only adequately qualified operators (e.g. machine operators, chemists, chemical technical assistants, physicist etc.) may work with the accessory. Operators must be trained before handling the accessory. Verify that the operators have read and understood the installation guide. Define precise responsibilities for the operators. Personal protective equipment must be provided to the operators.

- The responsible body must install a condensation water / thermal fluid drip tray below the temperature control unit (including accessory).
- The operator must check whether national regulations require the mandatory installation of a drain tray for the installation area of the temperature control unit (including accessory) / the entire system.

- Our temperature control unit (including accessory) complies with all applicable safety standards.
- Your system, which uses our temperature control unit (including accessory), must be as safe.
- The responsible body must design the system to ensure it is safe.
- Huber is not responsible for the safety of your system. The responsible body is responsible for the safety of the system.
- Whilst the temperature control unit (including accessory) provided by Huber meets all the applicable safety standards, integration into a system may give rise to hazards that are characteristic of the other system's design and beyond the control of Huber.
- It is the responsibility of the system integrator to ensure that the overall system, into which this temperature control unit (including accessory) is integrated, is safe.
- The **>Mains isolator<** [36] on the temperature control unit (if present) may be provided with a facility to lock the main isolator in the off position to facilitate safe system installation and maintenance of the temperature control unit (including accessory). Accessories with their own power supply must **additionally** be disconnected from the mains power supply! It is the responsibility of the responsible body to develop any lock-out/tag-out procedure in accordance with local regulations (e.g. CFR 1910.147 for the US).
- The responsible body must ensure that the cooling water supply and drain pipes are laid frost-protected. The cooling water temperature must not fall below 3 °C. The cooling water supply must be heated if the ambient temperature is below 3 °C.

1.2.1.1 Proper disposal of resources and consumables

Do comply with all national disposal regulations applicable for you. Contact your local waste management company for any questions concerning disposal.

Material / Aids	Disposal / Cleaning
Packing material	Keep the packaging material for future use (e.g. transport).
Thermal fluid	Please refer to the safety data sheet of the thermal fluid used for information on its proper disposal. Use the original thermal fluid container when disposing it.
Filling accessories, e.g. beaker	Clean the filling accessories for reuse. Make sure that the materials and cleaning agents used are properly disposed of.
Aids such as towels, cleaning cloths	Tools used to take up spilled thermal fluid must be disposed of in the same fashion as the thermal fluid itself. Tools used for cleaning must be disposed of depending on the cleaning agent used.
Cleaning agents such as stainless steel cleaning agents, sensitive-fabrics detergents	Please refer to the safety data sheet of the cleaning agent used for information on its proper disposal. Use the original containers when disposing of large quantities of cleaning agents.
Consumables such as air filter mats, temperature control hoses	Please refer to the safety data sheet of the consumables used for information on their proper disposal.

1.2.2 Requirements for operators

Work on the temperature control unit / accessory is reserved for appropriately qualified specialists, who have been assigned and trained by the responsible body to do so. Operators must be at least 18 years old. Persons under the age of 18 years may operate the temperature control unit / accessory only under the supervision of a qualified specialist. The operator is responsible for other people within the unit's working range.

1.2.3 Obligations of the operators

Carefully read the installation guide before handling the temperature control unit / accessory. Always observe the safety instructions. Wear appropriate personal protective equipment (e.g. safety goggles, protective gloves, non-slip shoes) when operating the temperature control unit / accessory.

1.3 General information

1.3.1 Description of workstation

The workstation is located at the control panel in front of the temperature control unit. The workstation is determined by the customer's connected peripheries. Accordingly, it must be designed safe by the responsible body. The workstation design also depends on the applicable requirements of the German occupational health and safety regulations [BetrSichV] and the risk analysis for the workstation.

2 Commissioning

2.1 Unpacking

PROCEDURE

- Check for damage to the packaging. Damage can indicate property damage to the accessory.
- Check for any transport damage when unpacking the accessory.
- Always contact your forwarding agent regarding the settlement of claims.
- Follow the instructions on page 12, section »Proper disposal of resources and consumables« for the disposal of packaging material.

2.2 Ambient conditions

**Unsuitable ambient conditions/unsuitable installation
SERIOUS INJURY DUE TO CRUSHING**

- Comply with the requirements under sections »Ambient conditions« and »Installation conditions«.

INFORMATION

Make sure there is adequate fresh air available at the site for the accessory. The warm exhaust air must be able to escape upwards unhindered.

Use of the accessory is permitted only under normal ambient conditions in accordance with DIN EN 61010-1:2011:

- Use only indoors. The illuminance must be at least 300 lx.
- Installation elevation up to 2000 meters above sea level.
- Maintain wall and ceiling clearance for adequate air exchange (dissipation of waste heat, supply of fresh air for the accessory and work area). Do not operate the accessory from within the cardboard box or with an inadequately dimensioned bath as this inhibits the air exchange.
- Ambient temperature values are provided on the technical data sheet; to ensure trouble-free operation, compliance with the ambient conditions is mandatory.
- Relative humidity max. 80% up to 32 °C and up to 40 °C decreasing linearly to 50%.
- Short distance to supply connections.
- The accessory must not be installed so as to hinder or even prevent access to the disconnecting device (to the power supply).

Wall clearances

Side	Distance in cm
Top	free standing
Front	min. 10
Right	min. 10
Left	min. 10
Rear	min. 10

2.3 Installation conditions

- Allow the temperature control unit to acclimate for about 2 hours when changing from a cold to a warm environment (or vice versa). Do not turn on the temperature control unit before!
- Install upright, stable and without tilt.
- Use a non-combustible, sealed foundation.
- Keep environment clean: Prevent slip and trip hazards.
- Spilled/leaked thermal fluid must be discarded immediately and correctly. Follow the instructions for the disposal of thermal fluid and material on page 12 in section »Proper disposal of resources and consumables«.
- Observe the ambient conditions.

2.4 Recommended temperature control and cooling water hoses

CAUTION

Use of unsuitable/defective hoses and/or hose connections

INJURIES

- **Thermal fluid**
- Use appropriate hoses and/or hose connections.
- Check periodically for leaks and the quality of the hose and hose connections and take suitable measures (replace) as required.
- Isolate and protect temperature control hoses against contact/mechanical load.
- **Cooling water**
- Reinforced hoses must be used to satisfy tougher safety requirements.
- Shut off the cooling water supply to the temperature control unit even for shorter downtimes (e.g. overnight).

CAUTION

Hot or cold thermal fluid and surfaces

BURNS TO LIMBS

- Avoid direct contact with the thermal fluids or the surfaces.
- Wear your personnel protective equipment (e.g. temperature-resistant safety gloves, safety goggles, safety footwear).

CAUTION

Uncontrolled formation of ice at the connections and hoses of the thermal fluid circuit

SLIP AND TRIP HAZARD

- If the temperature is controlled in the minus range, ice forms at the hoses and connections of the thermal fluid circuit. This occurs by condensing and freezing of atmospheric humidity.
- Check the strength of the ice formation. If too much ice is formed, this increases the risk of the temperature control unit tipping over. Secure the temperature control unit against tipping if this is the case.
- Check the ground below the ice formation for condensation water. Collect the condensation water with a suitable container or thoroughly remove it at regular intervals. You thus prevent the danger of slipping caused by condensation.

To connect applications, use only temperature control hoses that are compatible with the thermal fluid used. When selecting temperature control hoses, also pay attention to the temperature range in which the hoses are to be used.

- We recommend you use only temperature-insulated temperature control hoses with your temperature control unit. The responsible body is responsible for the insulation of connection valves.
- We **exclusively** recommend reinforced hoses for connecting to the cooling water supply. Cooling water and insulated temperature control hoses can be found in the Huber catalogue under Accessories.

2.5 Wrench sizes and torques

Observe the proper wrench sizes for the thermal fluid connection at the accessory. The following table lists the thermal fluid connections and the resulting wrench sizes, as well as the torque values. Always perform a leak test afterwards and re-tighten the connections if required. The values of the maximum torque (see table) must **not** be exceeded.

Overview
wrench sizes and
torques

Connection	Sleeve nut wrench size	Connector wrench size	Recommended torques in Nm	Maximum torques in Nm
M16x1	19	17	20	24
M24x1.5	27	27	47	56
M30x1.5	36	32	79	93
	36	36	79	93
M38x1.5	46	46	130	153
G-thread (flat-sealing)	Adapt the torque to the material of the flat seal used. First hand-tighten the temperature control hose. When using adapters, do not overtighten the G-thread on the pump connection when connecting a temperature control hose. When connecting a temperature control hose to the adapter piece, secure the G thread against overwinding.			

2.6 Temperature control units with cooler

WARNING

**Open electrical cables below the accessory if the temperature falls below the dew point
DEATH FROM ELECTRICAL SHOCK BY WATER ENTRY INTO THE ELECTRIC LINES.**

- A temperature below the dew point may result in condensation in/on the accessory when and at the cooling water connections. The condensation is caused by high humidity at the cooling water-bearing components. The condensation exits the accessory at the bottom.
- Electrical lines directly below the accessory must be protected against liquid ingress.

CAUTION

The accessory is not connected to the cooling water return line at a thermal fluid temperature ≥ 100 °C and are actively switched on

SEVERE BURNS/SCALDING DUE TO HOT STEAM ESCAPING FROM THE COOLING WATER RETURN LINE

- Do not switch on the accessory at a thermal fluid temperature ≥ 100 °C.
- Do not install a shut-off valve in the cooling water return line in the building!
- The accessory may only be switched on if: a) The accessory has been connected to the cooling water return line in the building. b) The thermal fluid temperature is < 100 °C.
- Hot steam with a relatively high pressure can escape from the return line of the cooling water connection according to the water pressure / water flow rate.
- Only use suitable and temperature-stable hoses (e.g. reinforced hoses) and screw connections.
- Use suitable measures (e.g. fixed pipes, fastening the hoses) to ensure that a water hammer in the return line does not pose a risk.
- Very hot water (95 °C $> T > 60$ °C) can still escape for some time after the vapor phase.
- The recommended water inlet temperature and the recommended water pressure can be found in the data sheet from page 28 in the section »Annex«.

CAUTION

Use of unsuitable/defective hoses and/or hose connections

INJURIES

- **Thermal fluid**
- Use appropriate hoses and/or hose connections.
- Check periodically for leaks and the quality of the hose and hose connections and take suitable measures (replace) as required.
- Isolate and protect temperature control hoses against contact/mechanical load.
- **Cooling water**
- Reinforced hoses must be used to satisfy tougher safety requirements.
- Shut off the cooling water supply to the temperature control unit even for shorter downtimes (e.g. overnight).

CAUTION

Extremely hot / cold surfaces, connections and thermal fluids

BURNS/FREEZING OF LIMBS

- Surfaces, connections and tempered thermal fluids can be extremely hot or cold depending on the operating mode.
- Avoid direct contact with surfaces, connections and thermal fluids!
- Wear your personnel protective equipment (e.g. temperature-resistant safety gloves, safety goggles).

NOTE

No protection against corrosion

DAMAGE TO THE TEMPERATURE CONTROL UNIT

- The addition of anti-corrosion agents is mandatory if salts (chlorides, bromide) have been added to the water circuit.
- Ensure that the materials used in the cooling water circuit are resistant with respect to the cooling water. See the data sheet from page 28 in section »Annex« for information on the materials used.
- Take suitable measures to maintain the warranty conditions.
- For information about water quality, see www.huber-online.com.

NOTE

Usage of un-filtered river/sea or ocean water as cooling water

DAMAGE TO THE TEMPERATURE CONTROL UNIT

- Un-filtered river or sea water is not suitable for use as cooling water due to its contaminants.
- Use drinking water or filtered river or sea water for cooling.
- The particle size of suspended solids must not exceed 5 microns.
- Sea water must not be used for water cooling.
- For information about water quality, see www.huber-online.com.


NOTE

Leaking cooling water connections

DAMAGE BY ROOM FLOODING


- Slowly open the building-side shut-off valves of the cooling water supply and return line.
- If water leaks from the cooling water connections: shut off the cooling water supply and return line immediately.
- Provide leakproof cooling water connections.

Overview

- [1] Circulation flow
- [2] Circulation return
- [13] Cooling water inlet
- [14] Cooling water outlet

Installation example

INFORMATION

The accessory scales at cooling water temperatures > 60 °CC according to the water quality. De-scale regularly and on time. Please note on page 12 in the section »**Proper disposal of resources and consumables**«.

When connecting to the public water network, observe all legal regulations and those specified by your water supplier.

The responsible body must ensure that the cooling water supply and discharge pipes are laid so they are protected from frost. The cooling water temperature must not fall below 3 °C. At ambient temperatures below 3 °C, the cooling water supply must be heated.

The minimum pressure differential in the cooling water circuit and the recommended cooling water inlet temperature can be found on the data sheet (from page 28 in the section »**Annex**«).

PROCEDURE

- Connect the >**Circulation flow**< [1] on the accessory with the >**Circulation return**< [2] on the temperature control unit.
- Connect the >**Cooling water inlet supply**< [13] on the accessory with the building's water supply line. A seal must be used for this.
- Connect the >**Cooling water outlet supply**< [14] on the accessory with the building's water return line. A seal must be used for this.
- Only use the accessory at a thermal fluid temperature < 100 °C.
- If thermal fluid temperatures are high (between 60 °C und 95 °C) ensure sufficient through flow in order to prevent excessive scaling of the accessory and the return line.
- Shut off the cooling water supply during heating. You will save energy and achieve shorter heating times.
- Check the connections for leaks.

2.7 Connecting an externally closed/open application

PROCEDURE

- Connect your external application (see operation manual of the temperature control unit used).

3 Function description

3.1 Function description of the accessory

3.1.1 General functions

The **accessory** can only be used in combination with a **Huber temperature control unit**. It is ideal for counter-cooling. The **accessory** is installed in the return line (>Circulation return< [2]) of the temperature control unit.

4 Setup mode

4.1 Filling and draining the accessory

PROCEDURE

- For filling, venting and draining the temperature control unit, proceed as described in the operation manual of the temperature control unit.
- When removing the accessory, watch out for any thermal fluid and cooling water residues it may contain.
- Collect the residues and dispose of them properly.

5 Normal operation

NOTE

When the temperature control unit is switched off, the thermal fluid temperature in the accessory is above/below room temperature

DAMAGE TO THE ACCESSORY

- Use the temperature control unit to temper the thermal fluid in the accessory to room temperature (20 °C).
- Do not close the shut-off valves in the thermal fluid circuit.

5.1 Counter-cooling

5.1.1 Start counter-cooling

PROCEDURE

- Observe all warnings, information and instructions on page 16 in the section »**Temperature control units with cooler**«!
- Open the building's shut-off valve in the cooling water inlet pipe.

5.1.2 End counter-cooling

PROCEDURE

- Use the temperature control unit to set the temperature of the thermal fluid to room temperature (20 °C).
- Stop the thermoregulation.
- Close the building's shut-off valve in the cooling water inlet pipe.

6 Service/maintenance

6.1 Maintenance

Cleaning/maintenance while the temperature control unit is operating

MORTAL DANGER FROM ELECTRIC SHOCK

- Stop an ongoing temperature control process.
- Turn off the temperature control unit.
- Also disconnect the temperature control unit from the power supply.

Carrying out maintenance work not described in this installation guide

DAMAGE TO THE ACCESSORY

- For maintenance work not described in the installation guide, contact Huber company.
- Maintenance work not described in this installation guide is reserved for qualified specialists trained by Huber.
- Only perform the following maintenance work on the accessory yourself.

6.1.1 Function check and visual inspection

Inspection intervals

Cooling*	Description	Maintenance interval	Comment	Person responsible
A/W	Visually inspect hoses and hose connections	Prior to switching on the temperature control unit	Exchange leaking hoses and hose connections prior to switching on the temperature control unit. Please see on page 22 in the section »Replacing temperature control or coolant hoses«.	Responsible body and/or operators
A/W	Thermal fluid inspection	As required	–	Responsible body and/or operators
A/W	Inspect the accessory for damage	Every 12 months or after a change of location	–	Responsible body and/or operators
W	Inspect cooling water quality	Every 12 months	Descale the cooling water circuit as required. Documentation on water quality available at: www.huber-online.com	Responsible body and/or operators

*A = Air cooling; W = Water cooling; U = Applicable only for Unistats

6.1.2 Replacing temperature control or coolant hoses

Replace defective temperature control and/or coolant hoses **before** turning on the temperature control unit.

6.1.2.1 Replacing temperature control hoses

PROCEDURE

- To drain the temperature control unit/accessory, proceed as described on page 20 in the section »Filling and draining the accessory«.
- Replace defective temperature control hoses. When disposing of them, observe page 12 in the section »Proper disposal of resources and consumables«.
- Connect your external application (see operation manual of the temperature control unit used).
- To fill and vent the temperature control unit/accessory, proceed as described on page 20 in the section »Filling and draining the accessory«.
- Restart the temperature control unit in normal mode.

6.1.2.2 Replacing coolant hoses

PROCEDURE

- Drain the cooling water as described on page 25 in section »Draining the cooling water«.
- Replace the defective coolant hoses. When disposing of them, observe on page 12 section »Proper disposal of resources and consumables«.
- Reconnect the temperature control unit to the building's cooling water supply. Proceed as described on page 16 in section »Temperature control units with cooler«.
- Restart the temperature control unit in normal mode.

6.2 Thermal fluid inspection, replacement and circuit cleaning

Check and change the thermal fluid and clean the thermal fluid circuit as described in the operation manual of the temperature control unit.

6.3 Cleaning the surfaces

Extremely hot / cold surfaces, connections and thermal fluids

BURNS/FREEZING OF LIMBS

- Surfaces, connections and tempered thermal fluids can be extremely hot or cold depending on the operating mode.
- Avoid direct contact with surfaces, connections and thermal fluids!
- Wear your personnel protective equipment (e.g. temperature-resistant safety gloves, safety goggles).

A standard stainless steel cleaning agent is suitable for cleaning the stainless steel surfaces. Carefully clean painted surfaces (damp only) using a solution of sensitive-fabrics detergent. Follow the instructions on page 12, section »Proper disposal of resources and consumables« for the disposal of cleaning agents and material.

6.4 Decontamination/repairs

Returning an accessory for repair that was not decontaminated

PHYSICAL INJURY AND PROPERTY DAMAGE CAUSED BY HAZARDOUS MATERIALS IN OR ON THE ACCESSORY

- Carry out appropriate decontamination.
- The decontamination process depends on the type and quantity of the materials used.
- Consult the relevant safety data sheet.
- You will find a prepared return receipt at www.huber-online.com.

As the responsible body you are responsible for carrying out decontamination **before** third-party personnel come into contact with the temperature control unit / accessory. Decontamination must be carried out **before** the temperature control unit / accessory is returned for repair or inspection. Attach a clearly visible written notice stating that the temperature control unit / accessory has been decontaminated.

To simplify the process, we have prepared a form for you. This is available for download at www.huber-online.com.

7 Shutting down

7.1 Safety instructions and basic principles

CAUTION

Non-compliance with the safety data sheet for the thermal fluid to be used

INJURIES

- Risk of injury to the eyes, skin, respiratory tract.
- The safety data sheet for the thermal fluid to be used must be read prior to using it and its content must be respected.
- Observe the local regulations/work instructions.
- Wear your personal protective equipment (e.g. temperature-resistant safety gloves, safety goggles, safety footwear).
- Danger of slipping because floor and work area are contaminated. Clean the work station and follow the instructions for the disposal of thermal fluid and material on page 12 in section »Proper disposal of resources and consumables«.

CAUTION

Hot or very cold thermal fluid

SEVERE BURNS/FROSTBITE OF LIMBS

- Before draining, ensure that the thermal fluid has room temperature (20 °C).
- If, at this temperature, the thermal fluid is too viscous to be drained: Control the temperature of the thermal fluid for a few minutes until the viscosity will allow drainage. Never control the temperature of the thermal fluid when the Drain is open.
- Danger of burns when draining thermal fluids at temperatures above 20 °C.
- Wear your Personal Protective Equipment for draining.
- Only drain using suitable drainage hose and collecting container. These must be compatible with the thermal fluid and its temperature.

CAUTION

Extremely hot / cold surfaces, connections and thermal fluids

BURNS/FREEZING OF LIMBS

- Surfaces, connections and tempered thermal fluids can be extremely hot or cold depending on the operating mode.
- Avoid direct contact with surfaces, connections and thermal fluids!
- Wear your personnel protective equipment (e.g. temperature-resistant safety gloves, safety goggles).

INFORMATION

All safety instructions are important and must be followed during working operations according to the installation guide!

7.2 Switch-off

PROCEDURE

- Switch off the temperature control unit as described in the operation manual of the temperature control unit.

7.3 Draining the cooling water

CAUTION

Pressurized cooling water connections

RISK OF INJURY

- Wear your personnel protective equipment (e.g. safety goggles).
- Carefully open the cooling water connection. Open slowly (1-2 signal edges) and drain the cooling water slowly.

NOTE

The building's shut-off valve is not closed

PHYSICAL DAMAGE DUE TO FLOODING OF THE ROOMS

- Close the building's shut-off valve in the cooling water inlet pipe.
- Place a collecting container below the >Cooling water inlet< [13] and >Cooling water outlet< [14].

PROCEDURE

- Close the building's shut-off valve in the cooling water inlet pipe.
- Place a collecting container below the >Cooling water inlet< [13] and >Cooling water outlet< [14].
- Open the >Cooling water inlet< [13] and the >Cooling water outlet< [14]. The cooling water will begin to drain out. Allow all the cooling water to drain out to prevent the risk of freezing during transport and storage!

7.4 Draining the accessory

CAUTION

Hot or very cold thermal fluid

SEVERE BURNS/FROSTBITE OF LIMBS

- Before draining, ensure that the thermal fluid has room temperature (20 °C).
- If, at this temperature, the thermal fluid is too viscous to be drained: Control the temperature of the thermal fluid for a few minutes until the viscosity will allow drainage. Never control the temperature of the thermal fluid when the Drain is open.
- Danger of burns when draining thermal fluids at temperatures above 20 °C.
- Wear your Personal Protective Equipment for draining.
- Only drain using suitable drainage hose and collecting container. These must be compatible with the thermal fluid and its temperature.

PROCEDURE

- To drain the accessory proceed as described on page 20 in the section »Filling and draining the accessory«.
- When removing the accessory, watch out for any thermal fluid and cooling water residues it may contain.
- Collect the residues and dispose of them properly.

7.5 Separating the accessory from the temperature control unit

PROCEDURE

- Please note: In order to disconnect the accessory from the temperature control unit, first discharge the cooling water and drain the accessory
- Remove the temperature control hose from the accessory's >Circulation return< [2].
- Disconnect the >Circulation flow< [1] on the accessory from the >Circulation return< [2] on the temperature control unit.

7.6 Packing

Use the original packaging wherever possible! Further information can be found on page 14 in section »Unpacking«.

7.7 Shipping

NOTE**Improper transport of accessory****PROPERTY DAMAGE**

➤ Comply with all requirements in this section to avoid damage to the accessory.

- Always use the original packaging for transport.
- Protect attachments from damage during transport!
- Secure with tensioning belts/lashing straps rated for the weight concerned.
- Additionally secure (depending on model) with plastic film, cardboard and straps.

7.8 Disposal

NOTE**Improper disposal****ENVIRONMENTAL DAMAGE**

➤ Spilled or leaked thermal fluid must be disposed of immediately and properly. Follow the instructions for the disposal of thermal fluid and auxiliary material on page 12 in section »Proper disposal of resources and consumables«.

Huber temperature control units and Huber accessories are made of high quality, recyclable materials. For example: Stainless steel 1.4301 / 1.4401 (V2A), copper, nickel, FKM, Perbunan, NBR, ceramic, carbon, Al-Oxid, red brass, brass, nickel-plated brass and silver solder. Proper recycling of the temperature control unit and accessories can actively help reduce CO₂ emissions in the production of these materials. Follow the laws and regulations of your jurisdiction when disposing material.

7.9 Contact data

INFORMATION

Please contact your supplier and/or local dealer **before** you return your accessories. The contact information can be found "Contact" on our home page www.huber-online.com. Please have the serial number of the accessories ready. The serial number can be found on the rating plate of the accessories.

7.9.1 Telephone number: Customer Support

If your country is not mentioned in the list below: The responsible service partner can be found on our homepage www.huber-online.com under the heading „Contact“.

- Huber Deutschland: +49 781 9603 244
- Huber China: +86 (20) 89001381
- Huber India: +91 80 2364 7966
- Huber Ireland: +44 1773 82 3369
- Huber Italia: +39 0331 181493
- Huber Swiss: +41 (0) 41 854 10 10
- Huber UK: +44 1773 82 3369
- Huber USA: +1 800 726 4877 | +1 919 674 4266

7.9.2 Telephone number: Sales

Telephone: +49-781-9603-123

7.9.3 Email address: Customer Support

Email: support@huber-online.com

7.10 Certificate of Compliance

Please read page 23, section »Decontamination/repairs«.

8 Annex

8.1 Technical data

Overview	Designation	Description
	Permissible temperature range	See rating plate on the accessory
	Max. thermal fluid temperature when switching on the accessory	See rating plate on the accessory
	Cooling water quality	On page 16 in the section »Temperature control units with cooler«
	Cooling water hoses	On page 15 in the section »Recommended temperature control and cooling water hoses«
	Max. cooling water pressure	See rating plate on the accessory
	Min. cooling water differential pressure	See rating plate on the accessory
	Cooling water inlet temperature	3 °C ... 40 °C
	Thermal fluid connections	M24x1.5
	Cooling water connections	M16x1
	Material	Stainless steel 1.4301

Inspired by **temperature** designed for you

Peter Huber Kältemaschinenbau SE
Werner-von-Siemens-Str. 1
77656 Offenburg / Germany

Telefon +49 (0)781 9603-0
Telefax +49 (0)781 57211

info@huber-online.com
www.huber-online.com

Technischer Service: +49 (0)781 9603-244

-125 °C ... +425 °C

huber