

# Unistat 405w

**Unistat 405w controls De Dietrich Process Systems 16l Glass Lined Stainless Steel reactor**



**Requirement**

This case study demonstrates the minimum achievable process temperature, process temperature control and process temperature stability when a Unistat 405w is connected with a De Dietrich Process Systems 16-liter Glass Lined Stainless Steel reactor.

**Method**

The Unistat 405w was connected to the 16l DDPS QVF GLSS reactor with 1 x 1-meter vacuum insulated hose and 1 x 1.5-meter vacuum insulated hose. The process mass was simulated with 12l of Huber's "DW-Therm" inside the reactor. Under "Process control" from a Pt100 located inside the process mass, the temperature of the process was cycled through various set-points and the results recorded using Huber's "Spy Service" software via a USB thumb drive inserted in the USB interface on the Pilot ONE controller.

**Setup details**


- Temperature range: -45°C...+250°C
- Heating power: 3.0 kW
- Hoses: 2 x M30 Metal Insul.
- HTF: M90.055/170.03
- Reactor: DDPS GLSS 16L50 liter
- Reactor content: 12 l DW-Therm
- Control: process
- Stirrer speed: 200 rpm
- Amb. temperature: +22°C

## Results

### 1. Performance


The graphic below shows the speed and accuracy as the Unistat 405w cycles the process temperature between -20°C and 100°C.

Start (°C)	End (°C)	Time Taken	Ramp Rate
20°C	-20°C	88 Minutes	0.45 K/Min
100°C	-20°C	160 Minutes	0.75 K/Min
-20°C	20°C	38 Minutes	1 K/Min
20°C	100°C	100 Minutes	0.8 K/Min
100°C	20°C	130 Minutes	0.6 K/Min



## 2. Minimum Achievable Process Temperature

The graphic shows that the minimum achievable process temperature was -28.9°C.



## 3. Stability

The graphic below demonstrates the Unistat 405w's ability to hold the process at 20°C +/- 0.01K over a period of 12-hours.

