

Unistat® 405w

Cooling a DDPS 2-litre glass reactor to T_{min}

Requirement

This case looks at the minimum temperature that a Unistat 405w can take the process in a 2-litre DDPS jacketed reactor under "process" control.

Method

The Unistat 405w was connected to the reactor using two 1-metre insulated metal hoses. The reactor was filled with 1.5 litre of "M90.055.03", a silicon based HTF.

Results

Process temperature reaches $-20\text{ }^{\circ}\text{C}$ from $20\text{ }^{\circ}\text{C}$ (40 K) within 23 minutes (1.7 K/min.) and asymptotes at $-34\text{ }^{\circ}\text{C}$ after 1 hour.

Setup details

Unistat® 405w & DDPS reactor

Temperature range: $-45\text{...}250\text{ }^{\circ}\text{C}$
 Cooling power: $1.3\text{ kW @ }250\text{...}0\text{ }^{\circ}\text{C}$
 $0.7\text{ kW @ }-20\text{ }^{\circ}\text{C}$
 Heating power: $1.5\text{ kW / }3\text{ kW}$
 Pump speed: 3300 rpm
 Hoses: $2 \times 1\text{ m; M24} \times 1.5$
 (#9325)
 HTF: DW-Therm (#6479)
 Reactor: 2-litre jacketed glass reactor
 Reactor contents: $1.5\text{ litre M90.055.03}$
 (#6259)
 Reactor stirrer speed: 115 rpm
 Control: process

